

INCUBATORS

- Dual Temperatures (37°C or 57°C)
- Pulsing light indicates temperature is maintained
- Illuminated power on/off
- Highly reliable dry heat block for even heat distribution
- Thermometer NIST traceable



Incubators are available as follows:

	<u>Cat. No.</u>
14 Cavity Incubator	I/15
Power Pack US	14/P1*
Power Pack European	14/P2*
Power Pack Australian	14/P3*
Thermometer	14/T1*

*Power Packs and thermometers are included with the incubator.

Rev. 2
Date: 3/26/04
SGM Part #: 7814