



SPORE SUSPENSIONS

7953* *G. stearothermophilus*; 9372* *B. atrophaeus*;
27142* *B. pumilus*; 35021* *B. subtilis* "5230"

- Known concentration of spores suspended in 40% ethanol or aqueous solution
- Each vial contains 10 mL of suspension
- Concentrations are specified per 0.1 mL

Certification:

Spore suspensions are certified for population, D-value, and Z-value (steam). Spore suspensions are manufactured in compliance with USP and ISO 11138.

Spore suspensions are available as follows:

Spore suspension	Cat. No.
<i>G. stearothermophilus</i> 10 ³	SSS/3
<i>G. stearothermophilus</i> 10 ⁴	SSS/4
<i>G. stearothermophilus</i> 10 ⁵	SSS/5
<i>G. stearothermophilus</i> 10 ⁶	SSS/6
<i>B. atrophaeus</i> 10 ⁴	SSG/4
<i>B. atrophaeus</i> 10 ⁵	SSG/5
<i>B. atrophaeus</i> 10 ⁶	SSG/6
<i>B. atrophaeus</i> 10 ⁷	SSG/7
<i>B. pumilus</i> 10 ⁴	SSR/4
<i>B. pumilus</i> 10 ⁵	SSR/5
<i>B. pumilus</i> 10 ⁶	SSR/6
<i>B. pumilus</i> 10 ⁷	SSR/7
<i>B. subtilis</i> "5230" 10 ⁶	SS5230/6
<i>B. subtilis</i> "5230" 10 ⁷	SS5230/7

Other organisms available upon request.

Storage: 40% Ethanol: Freezer (between -25° and -10°C)
Aqueous: Refrigerate (2-8°C)

Expiration: 24 months from date of manufacture (stored in ethanol)
12 months from date of manufacture (stored in aqueous)

* Culture is traceable to a recognized culture collection identified in USP and ISO 11138.